

EMS WHOLESALE WAREHOUSE

2014

*Our mission is customer satisfaction through innovative and cost-effective products that improve health care delivery for caregivers and those they serve. We are proud of our current line of **MEDSOURCE** INTERNATIONAL CO.,LTD. brand products and we pride ourselves on the unique ability to offer new items that meet your individual needs as a customer.*

***MEDSOURCE** INTERNATIONAL CO.,LTD., is a international sales company in SHANGHAI, CHINA. also, we do have our professional sourcing company (**Medsource** China) over 11 years. and we (**Medsource** Factory, Inc.) do manufacturer lot's of medical products by ourself. In last 10 more years, we expand our products to 20 more countries, such as USA, Canada, UK, Ireland, Germany, France, Poland, Italy, Portugal, Australia, Malaysia, New zealand, Singapore, Brazil, Chile, Peru, Mexico, Puerto Rico ...*

MEDSOURCE INTERNATIONAL



>>> INFANT CPR TRAINING MANIKIN

This simulator demonstrates palpation, embryo heart sound auscultation, pelvic measurement and mammary nursing training. Uterus can be inflated to lifelike through air filling. The speed of fetal heart sound and volume can be adjusted.



ITEM #	DESCRIPTION	QTY	MOQ	FOB PRICE
MS-TM100	Infant CPR Training Manikin, PVC	Each	5	\$671.90

>>> HALF-BODY CPR MANIKIN

Half-body CPR Training Manikin

Features :

1. Simulating standard open airway
2. External breast compression: display device and alarm device
 - a. indicator light display of correct and wrong compression ; alarm of wrong compression
 - b. intensity display of correct (4-5 cm) and wrong (<4cm or >5cm) compression; alarm of wrong compression.
3. Artificial respiration (inhalation) : display device and alarm device
 - a. Inhalation <500-600ml or > 600ml ,wrong indicator light display and alarm prompting ;Inhalation between 500-600ml right indicator light display
 - b. Indicator light display open airway
 - c. Inhalation too quickly or too much result in air entering into stomach ; indicator light display and alarm prompting.
- 4 Ratio of compression and artificial respiration: 30:2 (one or two person).
5. Operating cycle: one cycle includes five times of 30:2 ratio of compression and artificial respiration.
6. Operation frequency: 100 times per minute (the latest international standard)
7. Operation methods: exercise operation
8. Examination of pupil response: mydriasis and myosis
9. Examination of carotid response : pinch the pressure ball by hand and simulate carotid pulse
10. Working Conditions: output power is 6V transformed from 220V external power through manostat.



ITEM #	DESCRIPTION	QTY	MOQ	FOB PRICE
MS-TM200	Half-body CPR Manikin	Each	5	\$278.90